


Aesculap Neurosurgery Diamond Knife




Incision with brilliant perfection

Aesculap Diamond Knife



Diamond Knife

Due to the properties of their unique blade material, diamond knives offer excellent cutting characteristics for highly clean, precise and force-free incision especially in neurosurgical applications.



REM view of a common scalpel blade



REM view of a diamond knife blade

- Blade made of natural diamond
- Superior mechanical stability & elasticity of the blade
- Sustained sharpness
- Excellently clean, precise and force-free incisions
- Protection mechanism for safe storage of the blade inside the handle
- Color coded TITANIUM handles
- Four different cutting geometries: round, retrograde, wedge or lancet blade



FD113D

Round blade, gold-colored, 7 facets
Length 205 mm



FD114D

Retro blade, copper-colored, 60°
Length 205 mm



FD115D

Wedge blade, black, 45°
Length 205 mm



FD116D

Lancet blade, bronze-colored, 60°
Length 205 mm



AESCULAP®

B | BRAUN
SHARING EXPERTISE

Aesculap AG & Co. KG

Am Aesculap-Platz
78532 Tuttlingen
Germany

Phone +49 7461 95-0
Fax +49 7461 95-2600

www.aesculap.de

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogues and pricelists and to take legal actions.